CE
0499

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PROFESSIONAL USE ONLY

Composition:

VITRO MOLAR Powder contains: Barium Aluminum Silicate, Dehydrated Polyacrylic Acid and Iron Oxide.

VITRO MOLAR Liquid contains: Polyacrylic Acid, Tartaric Acid and Distilled Water.

VITRO CONDICIONADOR contains: Polyacrylic Acid, Glycerin, Aerosil 200, Methylene Blue CI 52015 and Deionized Water.

VITRO MOLAR is a condensable, chemically activated glass ionomer cement, used in restorations of posterior teeth; it is radiopaque and offers good esthetic properties, easy manipulation, fast curing and it can be easily sculpted. Its jellification occurs through an acid/base reaction. In addition to the features that are common to other glass ionomers, i. e., biocompatibility, high chemical adhesion, low solubility, adhesion to wet surfaces and continuous fluoride release, **VITRO MOLAR** is also highly resistant to compression and bending, and has a hard, low-abrasive surface and high viscosity. It is ideal to be used in the atraumatic restorative treatment (ART) technique.

Indications:

VITRO MOLAR is recommended for final Class I, II, III and V restorations of deciduous teeth, permanent restorations or long-term restoration in non-impact Class I and II areas, and for intermediary restorations. It is recommended for atraumatic restorative treatment (ART) technique, for cavities involving one or more surfaces, restorations of deciduous/permanent teeth and filling cores. It can be used as an intermediary material for cavities under heavy stress (Class I and II) and in adequacy of oral environment in patients with or without activity.

Presentations:

Box containing: 1 bottle (5g) of **VITRO MOLAR** powder + 1 bottle (4mL) of **VITRO MOLAR** liquid + 1 measuring spoon. Boxes are available in the following colors: A1, A3 or C2.

Box containing: 1 bottle (10g) of **VITRO MOLAR** powder + 1 bottle (8mL) of **VITRO MOLAR** liquid + 1 measuring spoon + 1 mixing pad for spatulation. Boxes are available in the following colors: A1, A3 or C2.

Box containing: 1 bottle (10g) of **VITRO MOLAR** powder + 1 bottle (8mL) of **VITRO MOLAR** liquid + 1 bottle (10mL) of **VITRO CONDICIONADOR** + 1 measuring spoon + 1 mixing pad for spatulation. Boxes are available in the following colors: A1, A3 or C2.

Box containing: 1 bottle (12.5g) of **VITRO MOLAR** powder + 1 bottle (8.5mL) of **VITRO MOLAR** liquid + 1 measuring spoon + 1 mixing pad for spatulation. Boxes are available in the following colors: A1, A3 or C2.

Box containing: 1 bottle (10g) of **VITRO MOLAR** powder + 1 measuring spoon. Boxes are available in the following colors: A1, A3 or C2.

Box containing: 1 bottle (8g) of **VITRO MOLAR** liquid.

Box containing: 1 bottle (10mL) of **VITRO CONDICIONADOR**.

Storage and Conservation:

VITRO MOLAR must be stored between 15°C and 30°C. Do not expose to excessive heat. Keep the containers of powder and liquid tightly closed after use. Under these conditions, the product is valid for up to 3 years as of the manufacturing date. Do not use this product if shelf life had expired.

Instructions for Use:

The use of gloves is recommended during the procedure.

Attention: The patient must be informed on the risk to occur sensitivity in the place of the procedure. This sensitization can be minimized through the use of proper isolation (rubber sheet) and pulp protection.

Do not apply **VITRO MOLAR** on the pulp, or over a dentin layer less than 0.5 mm thick. In such cases, if endodontic treatment is not necessary, a calcium hydroxide cement should be previously applied instead.

Preparation

A healthy dentin can be completely preserved during the preparation. Because of the adherence of **VITRO MOLAR** onto the enamel and dentin, it is not necessary to perform any retention procedures. However, care must be taken to avoid the formation of dentin layers less than 0,5mm thick. After preparation, the cavity surface must be well cleaned with water and it must be dried without causing dehydration.

Tooth Conditioning

Before applying **VITRO MOLAR**, you should apply **VITRO CONDICIONADOR** over the cavity surface (enamel or dentin) for 10 seconds to remove the dentin mud. After that, rinse it thoroughly with water and carefully dry it with sterilized absorbent paper. Keep the surface wet. Note: Glass ionomers, unlike other materials, do not work well on a dehydrated surface.

Dosage

Gently shake the bottle, to release the powder. Fill the measuring spoon with powder and level it in the stopper. The proportion used should be of 1 measuring spoon of the powder to 1 drop of the liquid. Liquid should be dispensed holding the bottle in a vertical position. Any crystals present in the neck should be removed, to assure there are not air bubbles in the drop dispensed. Keep the powder and liquid bottles tightly closed after use.

Manipulation

With a metal or plastic spatula, prepare the mixture over an appropriate mixing pad or glass slab during 20 seconds. The ideal procedure is to add the powder into the liquid in 2 portions, adding the second portion just after complete homogenization of the first portion. Mix the powder with the liquid until a homogenous consistence is obtained.

Application

Apply **VITRO MOLAR** with **CENTRIX** syringe and Accudose Posterior HV tips. During application, do not remove the tip from the material already applied in order to prevent the formation of air bubbles. The appropriate technique for each kind of indication must be followed. Whenever possible, apply some pressure using an acetate matrix and your finger over the applied material until complete setting is obtained. Remove the matrix and immediately apply a varnish or **ALPHA BOND** on top of the restoration.

Finishing and Polishing

Wait 6 minutes before you start the finishing procedure. Depending on the area, the occlusal adjust and the excess material can be removed with finishing diamond burs for composite resins. Use low rotation. Polishing can be performed with silicone tips in low rotation. Apply a new varnish layer or **ALPHA BOND** and instruct the patient not to chew for about one hour.

Properties of Use

Times given below are valid for the specific temperature and humidity conditions described. Higher temperatures will **reduce** and lower temperatures will **increase** working time.

Mixing Time (21-23°C, 50-52% RH)	30 sec
Working Time (21-23°C, 50-52% RH)	> 2 min 10 sec
Net Setting Time (36-38°C, 90% RH)	2 min 30sec - 4 min

The product tends to adhere to metal instruments; for this reason, it is recommended that these instruments are cleaned with cold water before the residue dries off.

Product wastes must not be reused. They must be disposed as the local law in force.

Precautions:

Avoid the use of **VITRO MOLAR** in patients with previous history of allergy to glass ionomer cements. In case of allergy, you must discontinue the use immediately and advise the patient to see a doctor. The liquid and the prepared cement can be corrosive. Do not allow contact with soft tissues or the skin. In case of contact, you should immediately remove the material with hydrophilic cotton soaked with water and rinse thoroughly with water. Avoid contact with the eyes; in case of accidental contact, wash the eyes immediately with water and see a doctor. Do not mix the powder or liquid which are components of this product with any other glass ionomer-based product. This product must be used only by dentists and for the intended use described in the Instructions for Use. Avoid contact with eugenol-based products, since they can delay the action of glass ionomers.

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DFL Indústria e Comércio S.A.