## NATURAL SHADE

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**NATURAL SHADE**  
Light-cured Micro hybrid Restorative Composite  
English

**PROFESSIONAL USE ONLY**

**NATURAL SHADE** is a light-cured micro hybrid restorative resin, radiopaque with excellent consistency.

**NATURAL SHADE** is manufactured with nanotechnology, under strict quality control, which provides superior esthetic and functional results to the professionals.

**Composition:**

**NATURAL SHADE** contains:

- Bisphenol A Diglycidylmethacrylate (Bis-GMA)
- Triethyleneglycol Dimethacrylate (TEGDMA)
- Aluminum Borosilicate
- Silica
- Camphoroquinone
- Butylated hydroxytoluene (BHT)
- Pigments

**Technical Data According to ISO 4049:**

Average particle size: .....0.7 µm  
Range of particle size: .....0.04 - 2.8 µm  
Percentage by volume of total inorganic filler: ..... 60 - 64 %  
Percentage by weight of total inorganic filler: ..... 75 - 79 %

**NATURAL SHADE** is indicated for:

Restorations:

- Classes I, II, III, IV and V
- Direct Veneers
- Closing of Diastemas
- Periodontal splinting

**Instructions for use:**

- Remove completely the restoration and/or the caried tissue. If necessary, bevel the enamel margins;
- Cleanse the tooth to be treated and the adjacent teeth as well;
- Select the shades, ideally in broad day light;
- Isolate the cavity;
- If necessary, proceed with pulp protection by using a calcium hydroxide cement and/or Glass Ionomer VITRO FIL LC;
- Acid etch using ALPHA ACID or ALPHA ETCH, by 30 seconds in enamel and 15 seconds in the dentin;
- Rinse with water spray and remove the excess moisture avoiding dehydrating the dentin;
- If necessary, apply a transparent matrix band by using interdental wedges;
- Apply a layer of dentin adhesive NATURAL BOND DE, rubbing by 5 seconds. Wait for 20 seconds and apply gently a jet of air by 5 to 10 seconds;
- Apply a second layer of NATURAL BOND DE. After waiting for 20 seconds, light cure the adhesive by 20 seconds;
- Dose a portion of **NATURAL SHADE** directly onto the mixing pad. Apply **NATURAL SHADE** in the cavity in increments up to 2mm;

- Light cure each increment for 20 seconds, finalizing with a light cure of 40 seconds;
- Do the finishing in three stages:
  1. Finishing with diamond burs;
  2. Polishing with NATURAL POLUX siliconized polishers;
  3. Application of the glazer or surface sealant NATURAL GLAZE.

**Presentations:**

- Envelope containing 1 syringe with 4g of **NATURAL SHADE** + 1 bottle with 4mL of dentin adhesive **NATURAL BOND DE**.
- Kit containing 4 syringes with 4g of **NATURAL SHADE** + 1 bottle with 4mL of dentin adhesive **NATURAL BOND DE** + 1 syringe with 3g of acid conditioner **ALPHA ETCH**.
- Kit containing 5 syringes with 4g of **NATURAL SHADE** + 1 bottle with 4mL of dentin adhesive **NATURAL BOND DE** + 1 syringe with 3g of acid conditioner **ALPHA ETCH**.
- Kit containing 7 syringes with 4g of **NATURAL SHADE** + 1 bottle with 4mL of dentin adhesive **NATURAL BOND DE** + 1 syringe with 3g of acid conditioner **ALPHA ETCH**.
- Refil containing 1 syringe with 4g of **NATURAL SHADE**.
- Refil containing 2 syringes with 4g of **NATURAL SHADE**.
- Refil containing 3 syringes with 4g of **NATURAL SHADE**.
- Refil containing 4 syringes with 4g of **NATURAL SHADE**.
- Refil containing 5 syringes with 4g of **NATURAL SHADE**.

**NATURAL SHADE** is available in the following shades:

- **Enamel** A1; A2; A3; A3,5; A4; B1; B2; B3; C2; C3; C4; D3;
- **Enamel Bleach** CL;
- **Dentin** A2; A3; A3,5; B2; C2; UNIVERSAL;
- **Incisal**.

**Shelf Life and Storage:**

**NATURAL SHADE** must be stored between 2 and 30°C. Exposure to direct sunlight or other heat sources should be avoided.

**NATURAL SHADE** fulfills the intended purpose until 3 years after manufacturing date.

Do not use this product if shelf life has expired.

Product wastes must not be reused. They must be disposed as the local law in force.

**Precautions:**

**NATURAL SHADE** is for professional use only.

Allergic reactions to resin components or the amine component are possible in patients with a special disposition.

Avoid contact with the eyes.

Keep out of the reach of children.

This product must be used only by dentists and for the intended use above described.

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DFL Indústria e Comércio S.A.